

Anti-Smox Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nwm0

Item Number: CM0015548



Introduction

SMOX Rabbit pAb is a polyclonal antibody targeting spermine oxidase, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat SMOX. Ideal for polyamine metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	SMOX Rabbit pAb
Remarks/Alias	PAO; SMO; PAO1; PAOH; PAO-1; PAOH1; C20orf16; SMOX
Species	Human
Gene ID	54498
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-544 of human SMOX (NP_787033.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NWM0
Protein Weight	62kDa
Shipping	Ice bag