


Anti-Smpd1 / Asm Polyclonal Antibody For Wb, If/Icc, Elisa - P17405

Item Number: CM0028996



Introduction

Rabbit polyclonal antibody against human SMPD1/ASM (acid sphingomyelinase), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studies of sphingolipid metabolism, apoptosis, and Niemann-Pick disease.

[Learn More](#)

Parameter	Value
Product Name	SMPD1 / ASM Rabbit pAb
Remarks/Alias	ASM; NPD; ASMASE; SMPD1 / ASM
Species	Human
GeneID (Human)	6609
GeneID	6609
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-340 of human SMPD1 / ASM (NP_000534.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17405
Protein Weight	70kDa
Shipping	Ice bag