


Anti-Smpd1 (Ko Validated) Rabbit Monoclonal Antibody For Wb, If-F, Ihc-P, Elisa - P17405

Item Number: CM0007145



Introduction

KO-validated rabbit monoclonal antibody against human SMPD1 (acid sphingomyelinase). Suitable for WB, IF-F, IHC-P, and ELISA. Cross-reacts with human and rat. Ideal for sphingolipid metabolism and lysosomal storage disease research.

[Learn More](#)

Parameter	Value
Product Name	Anti-SMPD1 (KO Validated) Rabbit Monoclonal Antibody
Remarks/Alias	ASM; NPD; ASMASE; D1
Species	Human
Gene ID (Human)	6609
Gene ID	6609
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-340 of human SMPD1 (NP_000534.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P17405
Protein Weight	70 kDa
Shipping	Ice bag