

## Anti-Smpd2 Rabbit Polyclonal Antibody For Western Blot And Elisa - O60906

Item Number: CM0013373



### Introduction

SMPD2 Rabbit Polyclonal Antibody validated for WB and ELISA; detects human SMPD2 and cross-reacts with mouse and rat. Ideal for studying sphingomyelinase activity and lipid metabolism pathways.

[Learn More](#)

| Parameter        | Value                                                                                                    |
|------------------|----------------------------------------------------------------------------------------------------------|
| Product Name     | SMPD2 Rabbit pAb                                                                                         |
| Remarks/Alias    | ISC1; NSMASE; NSMASE1; SMPD2                                                                             |
| Species          | Human                                                                                                    |
| Gene ID (Human)  | 6610                                                                                                     |
| Gene ID          | 6610                                                                                                     |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 150-250 of human SMPD2 (NP_003071.2). |
| Source           | Rabbit                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                    |
| Application      | WB, ELISA                                                                                                |
| Cross Reactivity | Mouse, Rat                                                                                               |
| SWISS            | O60906                                                                                                   |
| Protein Weight   | 48kDa                                                                                                    |
| Shipping         | Ice bag                                                                                                  |