


Anti-Smpd3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ny59

Item Number: CM0014148



Introduction

Rabbit polyclonal antibody against human SMPD3 (nSMase2), suitable for WB and ELISA; cross-reacts with mouse and rat. Targets sphingomyelin phosphodiesterase 3 involved in ceramide-mediated apoptosis and bone mineralization.

[Learn More](#)

Parameter	Value
Product Name	SMPD3 Rabbit pAb
Remarks/Alias	NSMASE2; SMPD3
Species	Human
Gene ID	55512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 416-655 of human SMPD3 (NP_061137.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NY59
Protein Weight	71 kDa
Shipping	Ice bag