


Anti-Smpd3/Nsmase2 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Ny59

Item Number: CM0007714



Introduction

High-quality rabbit monoclonal antibody against human SMPD3/NSMase2 (SWISS: Q9NY59). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studies of sphingolipid metabolism, apoptosis, and bone mineralization.

[Learn More](#)

Parameter	Value
Product Name	SMPD3/NSMase2 Rabbit mAb
Remarks/Alias	NSMASE2; SMPD3/NSMase2
Species	Human
GeneID (human)	55512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 461-655 of human SMPD3/NSMase2 (NP_061137.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NY59
Protein Weight	71kDa
Shipping	Ice bag