


Anti-Smtn Polyclonal Antibody For Wb, Elisa - P53814

Item Number: CM0026317



Introduction

Smoothelin (SMTN) rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat SMTN, a cytoskeleton structural protein in smooth muscle. Ideal for vascular and contractile tissue studies.

[Learn More](#)

Parameter	Value
Product Name	SMTN Rabbit pAb
Remarks/Alias	SMTN
Species	Human
GeneID (Human)	6525
GeneID	6525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 65-325 of human SMTN (NP_001193946.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53814
Protein Weight	99kDa
Shipping	Ice bag