


Anti-Smu1 Polyclonal Antibody For Wb, Elisa - Q2Tay7

Item Number: CM0031089



Introduction

Rabbit polyclonal antibody targeting human SMU1, suitable for WB and ELISA. Detects SMU1 in human, mouse, and rat samples. Ideal for spliceosome, mitosis, and viral infection research.

[Learn More](#)

Parameter	Value
Product Name	SMU1 Rabbit pAb
Remarks/Alias	BWD; SMU-1; fSAP57; SMU1
Species	Human
Gene ID (Human)	55234
Gene ID	55234
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human SMU1 (NP_060695.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q2TAY7
Protein Weight	58kDa
Shipping	Ice bag