

Anti-Smug1 Polyclonal Antibody For Wb And Elisa - Q53Hv7

Item Number: CM0014660

Introduction

SMUG1 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and mouse SMUG1, a single-strand selective uracil DNA glycosylase involved in base excision repair. Ideal for DNA damage and repair studies.

[Learn More](#)

Parameter	Value
Product Name	SMUG1 Rabbit pAb
Remarks/Alias	FDG; UNG3; HMUDG; SMUG1
Species	Human
Gene ID (Human)	23583
Gene ID	23583
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human SMUG1 (NP_055126.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q53HV7
Protein Weight	30kDa
Shipping	Ice bag