


Anti-Smurf2 Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Hau4

Item Number: CM0006702



Introduction

Rabbit monoclonal antibody targeting human SMURF2 (Q9HAU4), reactive in human and mouse. Validated for WB, IHC-P, IP, ELISA. Designed with recombinant immunogen (aa 201-290). Ideal for studies on E3 ubiquitin ligase function and TGF-beta signaling.

[Learn More](#)

Parameter	Value
Product Name	SMURF2 Rabbit mAb
Remarks/Alias	SMURF2; E3 ubiquitin-protein ligase SMURF2
Species	Human
Gene ID (Human)	64750
Gene ID	64750
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 201-290 of human SMURF2 (Q9HAU4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HAU4
Protein Weight	86kDa
Shipping	Ice bag