

Anti-Smyd1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Nb12

Item Number: CM0020139

Introduction

Rabbit polyclonal antibody targeting SMYD1 (SET and MYND domain-containing protein 1), a histone H3K4 methyltransferase repressor critical for cardiomyocyte differentiation. Validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	SMYD1 Rabbit pAb
Remarks/Alias	BOP; KMT3D; ZMYND18; ZMYND22; SMYD1
Species	Human
GeneID (Human)	150572
GeneID	150572
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human SMYD1 (NP_938015.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NB12
Protein Weight	57kDa
Shipping	Ice bag