


Anti-Smyd2 Polyclonal Antibody For Wb/If/Elisa - Q9Nrg4

Item Number: CM0030636



Introduction

Rabbit polyclonal antibody targeting SMYD2 (KMT3C, HSKM-B) for WB, IF/ICC, ELISA. Cross-reacts with human and mouse. Ideal for studying histone methylation, p53 regulation, and cancer epigenetics.

[Learn More](#)

Parameter	Value
Product Name	SMYD2 Rabbit pAb
Remarks/Alias	KMT3C; HSKM-B; ZMYND14; SMYD2
Species	Human
GeneID (Human)	56950
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human SMYD2 (NP_064582.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NRG4
Protein Weight	50kDa
Shipping	Ice bag