


Anti-Smyd3 Monoclonal Antibody For Wb, If/Icc, Elisa - Q9H7B4

Item Number: CM0006979



Introduction

SMYD3 Rabbit Monoclonal Antibody for Western blot, immunofluorescence, and ELISA. Cross-reacts with human, mouse, rat SMYD3 protein (49 kDa). Suitable for chromatin modification and oncology research.

[Learn More](#)

Parameter	Value
Product Name	SMYD3 Rabbit mAb
Remarks/Alias	KMT3E; ZMYND1; ZNFN3A1; bA74P14.1; SMYD3
Species	Human
Gene ID (Human)	64754
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human SMYD3 (Q9H7B4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H7B4
Protein Weight	49kDa
Shipping	Ice bag