


Anti-Smyd3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H7B4

Item Number: CM0018187



Introduction

SMYD3 Rabbit Polyclonal Antibody validated for Western blot, IHC-P, IF/ICC, and ELISA in human, mouse, and rat. Used for histone methylation and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	SMYD3 Rabbit pAb
Remarks/Alias	KMT3E; ZMYND1; ZNFN3A1; bA74P14.1; SMYD3
Species	Human
Gene ID (Human)	64754
Gene ID	64754
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human SMYD3 (NP_001161212.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H7B4
Protein Weight	49kDa
Shipping	Ice bag