

Anti-Smyd5 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q6Gmv2

Item Number: CM0030464

Introduction

Rabbit polyclonal antibody targeting human SMYD5 (ZMYND23), a protein-lysine N-trimethyltransferase involved in translation elongation, histone modification, and embryonic stem cell regulation. Validated for immunofluorescence (IF/ICC) and ELISA.

[Learn More](#)

Parameter	Value
Product Name	SMYD5 Rabbit pAb
Remarks/Alias	RRG1; RAI15; NN8-4AG; ZMYND23; SMYD5
Species	Human
Gene ID (Human)	10322
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-400 of human SMYD5 (NP_006053.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q6GMV2
Protein Weight	47kDa
Shipping	Ice bag