

Anti-Snail Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95863

Item Number: CM0013823



Introduction

Rabbit polyclonal antibody against human Snail (SNAI1), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for epithelial-mesenchymal transition (EMT) and cancer research.

[Learn More](#)

Parameter	Description
Product Name	Snail Rabbit pAb
Remarks/Alias	SNA; SNAH; SNAIL; SLUGH2; SNAIL1; dJ710H13.1; Snail
Species	Human
GeneID (Human)	6615
GeneID	6615
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200-267 of human Snail (NP_005976.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95863
Protein Weight	29kDa
Shipping	Ice bag