

Anti-Snail Polyclonal Antibody (Snai1) For Wb, Ihc-P, If/Icc, Elisa - O95863

Item Number: CM0031065



Introduction

Rabbit polyclonal antibody targeting Snail/SNAI1, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studies on epithelial-mesenchymal transition and tumor progression.

[Learn More](#)

Parameter	Value
Product Name	Snail Rabbit pAb
Remarks/Alias	SNA; SNAH; SNAIL; SLUGH2; SNAI1; dj710H13.1; Snail
Species	Human
GeneID (Human)	6615
GeneID	6615
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 200-267 of human SNAI1 (NP_005976.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95863
Protein Weight	29kDa
Shipping	Ice bag