

Anti-Snap23 Polyclonal Antibody For Wb, Ihc-P, Elisa - O00161

Item Number: CM0018387



Introduction

High-quality rabbit polyclonal antibody against human SNAP23, validated for WB, IHC-P, and ELISA. Suitable for studying SNARE-mediated membrane fusion in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	SNAP23 Rabbit pAb
Remarks/Alias	SNAP-23; SNAP23A; SNAP23B; HsT17016; SNAP23
Species	Human
GeneID (Human)	8773
GeneID	8773
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human SNAP23 (NP_003816.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00161
Protein Weight	23kDa
Shipping	Ice bag