

Anti-Snap25 Polyclonal Antibody For Wb, If, Ihc, Elisa - P60880

Item Number: CM0016370



Introduction

High-quality rabbit polyclonal antibody against human SNAP25, validated for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, rat. Ideal for synaptic function and neurotransmitter release research.

[Learn More](#)

Parameter	Value
Product Name	SNAP25 Rabbit pAb
Remarks/Alias	SUP, RIC4, SEC9, SNAP, CMS18, RIC-4, SNAP-25, bA416N4.2, dJ1068F16.2, SNAP25
Species	Human
Gene ID (Human)	6616
Gene ID	6616
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-206 of human SNAP25 (NP_570824.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60880
Protein Weight	23kDa
Shipping	Ice bag