

Anti-Snap25 Rabbit Monoclonal Antibody For Wb, Elisa - P60880

Item Number: CM0006600



Introduction

High-quality rabbit monoclonal antibody against human SNAP25 (P60880), validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for neuroscience and synaptic function research.

[Learn More](#)

Parameter	Value
Product Name	SNAP25 Rabbit mAb
Remarks/Alias	SUP; RIC4; SEC9; SNAP; CMS18; RIC-4; SNAP-25; bA416N4.2; dJ1068F16.2; SNAP25
Species	Human
Gene ID (Human)	6616
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P60880
Protein Weight	23kDa
Shipping	Ice bag