

## Anti-Snap29 Rabbit Polyclonal Antibody For Wb, Elisa - O95721

Item Number: CM0031079



### Introduction

High-quality rabbit polyclonal antibody against human SNAP29 for WB and ELISA applications. Cross-reacts with mouse and rat. Ideal for research on autophagy, membrane fusion, and ciliogenesis. 29kDa.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | SNAP29 Rabbit pAb                                                                                                 |
| Remarks/Alias    | CEDNIK; SNAP-29; SNAP29                                                                                           |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 9342                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-258 of human SNAP29 (NP_004773.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | O95721                                                                                                            |
| Protein Weight   | 29kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |