

Anti-Snapin Polyclonal Antibody For Wb, Ihc-P, Elisa - O95295

Item Number: CM0014039

Introduction

Rabbit polyclonal antibody against human SNAPIN (SNAPAP, BLOC1S7), a SNARE-associated protein involved in lysosome biogenesis and vesicle trafficking. Validated for WB, IHC-P, ELISA with reactivity to human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	SNAPIN Rabbit pAb
Remarks/Alias	BLOS7; BORCS3; SNAPAP; BLOC1S7; SNAPIN
Species	Human
GeneID (Human)	23557
GeneID	23557
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-136 of human SNAPIN (NP_036569.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95295
Protein Weight	15kDa
Shipping	Ice bag