

Anti-Snd1 Rabbit Monoclonal Antibody [Ko Validated] For Wb, Ihc-P, If/Icc, Elisa - Q7Kzf4

Item Number: CM0005766



Introduction

Recombinant rabbit monoclonal anti-SND1 antibody [KO Validated] for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Targets the 102 kDa miRNA endonuclease and transcriptional coactivator p100/TDRD11. Ideal for gene regulation and virology research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SND1 Rabbit mAb
Remarks/Alias	p100; TDRD11; Tudor-SN
Species	Human
GeneID (Human)	27044
GeneID	27044
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-285 of human SND1 (NP_055205.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7KZF4
Protein Weight	102kDa
Shipping	Ice bag