

Anti-Snd1 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Ihc-P, If/Icc, Ip, Elisa - Q7Kzf4

Item Number: CM0030480



Introduction

Rabbit polyclonal antibody against human SND1, validated by knockout. Suitable for WB, IHC-P, IF/ICC, IP, ELISA. Recognizes human, mouse, rat SND1. Ideal for miRNA decay and transcriptional regulation research. High-quality IVD raw material from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SND1 Rabbit pAb
Remarks/Alias	p100; TDRD11; Tudor-SN; D1
Species	Human
GeneID (Human)	27044
GeneID	27044
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-285 of human SND1 (NP_055205.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7KZF4
Protein Weight	102kDa
Shipping	Ice bag