


Anti-Snip1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Tad8

Item Number: CM0018846



Introduction

Rabbit polyclonal antibody targeting human SNIP1, validated for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with rat. Essential for pre-mRNA splicing, NF-κB signaling, and miRNA biogenesis.

[Learn More](#)

Parameter	Specification
Product Name	SNIP1 Rabbit pAb
Remarks/Alias	PML1; PMRED; NEDHCS; SNIP1
Species	Human
Gene ID (Human)	79753
Gene ID	79753
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human SNIP1 (NP_078976.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8TAD8
Protein Weight	46kDa
Shipping	Ice bag