

Anti-Snrk Polyclonal Antibody For Wb And Elisa - Q9Nrh2

Item Number: CM0020933



Introduction

High-quality SNRK Rabbit Polyclonal Antibody validated for WB and ELISA, cross-reacts with human and mouse. Targets SNF-related kinase involved in hematopoietic proliferation and neuronal apoptosis. For research and IVD use.

[Learn More](#)

Parameter	Value
Product Name	SNRK Rabbit pAb
Remarks/Alias	HSNFRK; SNRK
Species	Human
GeneID (Human)	54861
GeneID	54861
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-240 of human SNRK (NP_060189.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NRH2
Protein Weight	84kDa
Shipping	Ice bag