

Anti-Snrnp200 Rabbit Polyclonal Antibody For Wb, Elisa - O75643

Item Number: CM0024755



Introduction

Rabbit polyclonal antibody against human SNRNP200 (U5-200KD helicase), validated for Western blot and ELISA. Suitable for spliceosome research and pre-mRNA splicing mechanism studies.

[Learn More](#)

Parameter	Value
Product Name	SNRNP200 Rabbit pAb
Remarks/Alias	BRR2, RP33, HELIC2, ASCC3L1, U5-200KD, SNRNP200
Species	Human
GeneID (Human)	23020
GeneID	23020
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1947-2136 of human SNRNP200 (NP_054733.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75643
Protein Weight	245kDa
Shipping	Ice bag