

Anti-Snrnp35 Polyclonal Antibody For Wb, Elisa - Q16560

Item Number: CM0018524



Introduction

Rabbit polyclonal antibody targeting human SNRNP35 (U11/U12 snRNP 35 kDa protein). Suitable for Western blot and ELISA, reacting with human and mouse samples. Ideal for spliceosome and mRNA splicing research.

[Learn More](#)

Parameter	Value
Product Name	SNRNP35 Rabbit pAb
Remarks/Alias	HM1; HM-1; U1SNRNPBP; SNRNP35
Species	Human
GeneID (Human)	11066
GeneID	11066
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human SNRNP35 (NP_073208.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16560
Protein Weight	29kDa
Shipping	Ice bag