


Anti-Snrnp40 Polyclonal Antibody For Wb, Elisa - Q96Di7

Item Number: CM0021370



Introduction

Highly specific rabbit polyclonal antibody against human SNRNP40 (40K, WDR57, PRP8BP), a core U5 snRNP spliceosome component. Validated for Western blot and ELISA, with broad mammalian cross-reactivity to mouse and rat. Ideal for splicing research.

[Learn More](#)

Parameter	Value
Product Name	SNRNP40 Rabbit pAb
Remarks/Alias	40K; SPF38; WDR57; PRP8BP; HPRP8BP; PRPF8BP; SNRNP40
Species	Human
GeneID (Human)	9410
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-357 of human SNRNP40 (NP_004805.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96DI7
Protein Weight	39kDa
Shipping	Ice bag