

Anti-Snrnp70 Polyclonal Antibody For Wb, Ihc-P, Elisa - P08621

Item Number: CM0025065



Introduction

Rabbit polyclonal antibody targeting human SNRNP70 (U1-70K). Validated for WB, IHC-P, ELISA; cross-reacts with Human, Mouse, Rat. Essential for spliceosome assembly and pre-mRNA splicing research.

[Learn More](#)

Parameter	Detail
Product Name	SNRNP70 Rabbit pAb
Remarks/Alias	RPU1, Snp1, U1AP, U170K, U1RNP, RNPU1Z, SNRP70, U1-70K, SNRNP70
Species	Human
GeneID (Human)	6625
GeneID	6625
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human SNRNP70 (NP_003080.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08621
Protein Weight	52kDa
Shipping	Ice bag