

Anti-Snrnp70 Rabbit Polyclonal Antibody For Wb, Elisa - P08621

Item Number: CM0015106



Introduction

Rabbit polyclonal antibody targeting human SNRNP70 (U1-70K), a key component of the U1 snRNP spliceosome. Validated for Western blot and ELISA in human, mouse, and rat samples. For IVD and research applications.

[Learn More](#)

Parameter	Value
Product Name	SNRNP70 Rabbit pAb
Remarks/Alias	RPUI1; Snp1; U1AP; U170K; U1RNP; RNPU1Z; SNRP70; U1-70K; SNRNP70
Species	Human
GeneID (Human)	6625
GeneID	6625
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human SNRNP70 (NP_003080.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08621
Protein Weight	52kDa
Shipping	Ice bag