

Anti-Snrpb2 Rabbit Polyclonal Antibody For Wb, Elisa - P08579

Item Number: CM0027801



Introduction

Rabbit polyclonal antibody targeting human SNRPB2 (U2B"), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for studying pre-mRNA splicing and spliceosome function. 25kDa.

[Learn More](#)

Parameter	Value
Product Name	SNRPB2 Rabbit pAb
Remarks/Alias	Msl1; U2B"; SNRPB2
Species	Human
GeneID (Human)	6629
GeneID	6629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human SNRPB2 (NP_003083.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08579
Protein Weight	25kDa
Shipping	Ice bag