


Anti-Snrpd1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62314

Item Number: CM0017212



Introduction

High-quality rabbit anti-SNRPD1 polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA applications. Reacts with human, mouse, and rat. Ideal for studying spliceosome assembly and pre-mRNA splicing mechanisms.

[Learn More](#)

Parameter	Value
Product Name	SNRPD1 Rabbit pAb
Remarks/Alias	SMD1; SNRPD; Sm-D1; HsT2456; SNRPD1
Species	Human
Gene ID (Human)	6632
Gene ID	6632
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-96 of human SNRPD1 (NP_008869.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62314
Protein Weight	13kDa
Shipping	Ice bag