

Anti-Snrpd3 Rabbit Polyclonal Antibody For Wb, Elisa - P62318

Item Number: CM0018272



Introduction

Rabbit polyclonal antibody targeting SNRPD3 (Sm-D3), a core spliceosomal protein. Validated for Western blot and ELISA, with reactivity across human, mouse, and rat. Suitable for pre-mRNA splicing and histone mRNA processing research.

[Learn More](#)

Parameter	Value
Product Name	SNRPD3 Rabbit pAb
Remarks/Alias	SMD3; Sm-D3; SNRPD3
Species	Human
GeneID (human)	6634
GeneID	6634
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 8-121 of human SNRPD3 (NP_001265585)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62318
Protein Weight	14kDa
Shipping	Ice bag