

Anti-Snrpe Polyclonal Antibody For Wb, Ihc-P, Elisa - P62304

Item Number: CM0030661



Introduction

High-quality rabbit polyclonal antibody targeting human SNRPE (Small nuclear ribonucleoprotein E). Validated for Western blot, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for pre-mRNA splicing and spliceosome studies.

[Learn More](#)

Parameter	Value
Product Name	SNRPE Rabbit pAb
Remarks/Alias	SME; Sm-E; HYPT11; snRNP-E; SNRPE
Species	Human
Gene ID	6635
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human SNRPE (NP_003085.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62304
Protein Weight	11kDa
Shipping	Ice bag