

Anti-Snrpf Polyclonal Antibody For Wb, Ihc-P, Elisa - P62306

Item Number: CM0015242



Introduction

Rabbit polyclonal antibody targeting human SNRPF (small nuclear ribonucleoprotein F), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. For pre-mRNA splicing and snRNP studies.

[Learn More](#)

Parameter	Value
Product Name	SNRPF Rabbit pAb
Remarks/Alias	SMF; Sm-F; snRNP-F; SNRPF
Species	Human
GeneID (human)	6636
GeneID	6636
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-86 of human SNRPF (NP_003086.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62306
Protein Weight	10kDa
Shipping	Ice bag