

Anti-Snw1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13573

Item Number: CM0014964



Introduction

Rabbit polyclonal antibody against human SNW1 (SKIP, NCoA-62) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Target involved in pre-mRNA splicing and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	SNW1 Rabbit pAb
Remarks/Alias	Bx42; SKIP; FUN20; Prp45; SKIIP; SKIP1; PRPF45; NCOA-62; SNW1
Species	Human
GeneID (Human)	22938
GeneID	22938
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-536 of human SNW1 (NP_036377.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13573
Protein Weight	61kDa
Shipping	Ice bag