


Anti-Snx1 Rabbit Monoclonal Antibody For Wb, Elisa - Q13596

Item Number: CM0009529



Introduction

Rabbit monoclonal antibody against human sorting nexin-1 (SNX1/VP55). Recombinant immunogen (aa 1-100), validated for WB and ELISA. Detects a 59kDa band. Essential for investigating retromer-dependent cargo sorting and endosomal trafficking pathways.

[Learn More](#)

Parameter	Value
Product Name	SNX1 Rabbit mAb
Remarks/Alias	VP55; HsT17379; SNX1
Species	Human
GeneID (Human)	6642
GeneID	6642
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human SNX1 (Q13596).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13596
Protein Weight	59kDa
Shipping	Ice bag