


Anti-Snx1 Rabbit Pab For Wb - Q13596

Item Number: CM0025990



Introduction

Rabbit polyclonal antibody targeting human SNX1 (Sorting nexin-1), validated for Western blot, immunofluorescence, and ELISA. Detects a 59 kDa band and cross-reacts with mouse and rat. Suitable for intracellular trafficking and retromer studies.

[Learn More](#)

Parameter	Description
Product Name	SNX1 Rabbit pAb
Remarks/Alias	VPS5; HsT17379; SNX1
Species	Human
Gene ID (Human)	6642
Gene ID	6642
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human SNX1 (NP_003090.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13596
Protein Weight	59kDa
Shipping	Ice bag