

Anti-Snx15 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nrs6

Item Number: CM0030098



Introduction

Rabbit polyclonal antibody against human SNX15, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Detects a 38kDa protein involved in intracellular trafficking. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	SNX15 Rabbit pAb
Remarks/Alias	HSAF001435; SNX15
Species	Human
GeneID (Human)	29907
GeneID	29907
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-342 of human SNX15 (NP_037438.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRS6
Protein Weight	38kDa
Shipping	Ice bag