


Anti-Snx3 Rabbit Polyclonal Antibody For Wb And Elisa - O60493

Item Number: CM0028191



Introduction

Rabbit polyclonal antibody targeting human SNX3 (Sorting nexin-3), validated for WB and ELISA with cross-reactivity to mouse and rat. Suitable for studies on endosomal trafficking, retromer complex function, and WNT signaling.

[Learn More](#)

Parameter	Value
Product Name	SNX3 Rabbit pAb
Remarks/Alias	SDP3; Grd19; MCOPS8; SNX3
Species	Human
GeneID (Human)	8724
GeneID	8724
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-162 of human SNX3 (NP_003786.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60493
Protein Weight	19kDa
Shipping	Ice bag