


Anti-Snx32 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q86Xe0

Item Number: CM0029633



Introduction

SNX32 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Involved in endosomal trafficking and neurite outgrowth; ideal for neuroscience and cell biology studies.

[Learn More](#)

Parameter	Value
Product Name	SNX32 Rabbit pAb
Remarks/Alias	SNX6B; SNX32
Species	Human
GeneID	254122
Immunogen	Recombinant fusion protein corresponding to amino acids 218-285 of human SNX32 (NP_689973.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86XE0
Protein Weight	46kDa
Shipping	Ice bag