


Anti-Snx5 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y5X3

Item Number: CM0013646



Introduction

Rabbit polyclonal antibody targeting human SNX5, validated for WB and ELISA. Cross-reacts with mouse and rat. Useful for studying intracellular trafficking and retromer function.

[Learn More](#)

Parameter	Value
Product Name	SNX5 Rabbit pAb
Remarks/Alias	SNX5
Species	Human
GeneID (Human)	27131
GeneID	27131
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 145-404 of human SNX5 (NP_055241.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5X3
Protein Weight	47kDa
Shipping	Ice bag