


Anti-Snx7 Polyclonal Antibody For Wb, Elisa - Q9Unh6

Item Number: CM0027391



Introduction

Rabbit polyclonal antibody targeting human SNX7, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for endocytosis and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	SNX7 Rabbit pAb
Remarks/Alias	SNX7
Species	Human
GeneID (Human)	51375
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 85-330 of human SNX7 (NP_057060.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UNH6
Protein Weight	45kDa
Shipping	Ice bag