

Anti-Snx9 Polyclonal Antibody For Wb And Elisa - Q9Y5X1

Item Number: CM0018563



Introduction

Rabbit polyclonal antibody targeting human SNX9, validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat SNX9, relevant for endocytosis and membrane trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	SNX9 Rabbit pAb
Remarks/Alias	SDP1; WISP; SH3PX1; SH3PXD3A; SNX9
Species	Human
Gene ID (Human)	51429
Gene ID	51429
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human SNX9 (NP_057308.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5X1
Protein Weight	67kDa
Shipping	Ice bag