

Anti-Soat1 Monoclonal Antibody For Wb, Elisa - P35610

Item Number: CM0006898



Introduction

Rabbit monoclonal antibody targeting human SOAT1 (sterol O-acyltransferase 1, ACAT-1), validated for Western blot and ELISA. Ideal for research on cholesterol metabolism, lipoprotein assembly, and dietary cholesterol absorption.

[Learn More](#)

Parameter	Value
Product Name	SOAT1 Rabbit mAb
Remarks/Alias	ACAT; SOAT; STAT; ACACT; ACAT1; ACAT-1; SOAT1
Species	Human
GeneID	6646
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-104 of human SOAT1. (NP_003092.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P35610
Protein Weight	65kDa
Shipping	Ice bag