


Anti-Soat1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P35610

Item Number: CM0029108



Introduction

SOAT1 (ACAT-1) rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Targets cholesterol metabolism; 65 kDa. Suitable for lipid research.

[Learn More](#)

Parameter	Value
Product Name	SOAT1 Rabbit pAb
Remarks/Alias	ACAT, SOAT, STAT, ACACT, ACAT1, ACAT-1, SOAT1
Species	Human
GeneID	6646
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human SOAT1 (NP_003092.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35610
Protein Weight	65 kDa
Shipping	Ice bag