

Anti-Socs1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15524

Item Number: CM0025412

Introduction

Rabbit polyclonal antibody against human SOCS1 (JAB/SSI-1), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying JAK/STAT signaling and cytokine regulation.

[Learn More](#)

Parameter	Value
Product Name	SOCS1 Rabbit pAb
Remarks/Alias	JAB; CIS1; SSI1; TIP3; CISH1; SSI-1; TIP-3; AISIMD; SOCS-1; SOCS1
Species	Human
GeneID	8651
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human SOCS1 (NP_003736.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15524
Protein Weight	24kDa
Shipping	Ice bag