

Anti-Socs2 Monoclonal Antibody For Wb, Ihc-P, Elisa - O14508

Item Number: CM0004013



Introduction

Anti-SOCS2 rabbit monoclonal antibody validated for WB, IHC-P, ELISA; reacts with human and rat. Recognizes suppressor of cytokine signaling 2, a regulator of growth hormone receptor ubiquitination and degradation.

[Learn More](#)

Parameter	Value
Product Name	SOCS2 Rabbit mAb
Remarks/Alias	CIS2; SSI2; Cish2; SSI-2; SOCS-2; STATI2; SOCS2
Species	Human
GeneID (human)	8835
GeneID	8835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human SOCS2 (NP_001257400.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	O14508
Protein Weight	22 kDa
Shipping	Ice bag