

## Anti-Socs4 Polyclonal Antibody For Wb And Elisa - Q8Wxh5

Item Number: CM0026450



### Introduction

High-quality rabbit polyclonal antibody targeting human SOCS4 (also known as SOCS7), validated in WB and ELISA. Cross-reacts with mouse, suitable for cytokine signaling and EGFR degradation studies.

[Learn More](#)

| Parameter        | Value                                                                                                   |
|------------------|---------------------------------------------------------------------------------------------------------|
| Product Name     | SOCS4 Rabbit pAb                                                                                        |
| Remarks/Alias    | SOCS7; SOCS4                                                                                            |
| Species          | Human                                                                                                   |
| GeneID (Human)   | 122809                                                                                                  |
| GeneID           | 122809                                                                                                  |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 340-440 of human SOCS4 (NP_543143.1) |
| Source           | Rabbit                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                   |
| Application      | WB, ELISA                                                                                               |
| Cross Reactivity | Human, Mouse                                                                                            |
| SWISS            | Q8WXH5                                                                                                  |
| Protein Weight   | 51kDa                                                                                                   |
| Shipping         | Ice bag                                                                                                 |